




US 20170224267A1

(19) **United States**(12) **Patent Application Publication**  
**STRADER et al.**(10) **Pub. No.: US 2017/0224267 A1**(43) **Pub. Date: Aug. 10, 2017**(54) **PRICK TEST SINGLE USE STERILE VIAL  
AND METHOD***B65D 45/16* (2006.01)*B65B 3/00* (2006.01)*B65D 25/10* (2006.01)*B65D 25/20* (2006.01)(71) Applicant: **ROCA MEDICAL LTD.**, London (GB)(72) Inventors: **JAMES STRADER**, AUSTIN, TX  
(US); **JOVAN HUTTON PULITZER**,  
FRISCO, TX (US)(52) **U.S. Cl.**CPC ..... *A61B 5/411* (2013.01); *B65D 25/10*  
(2013.01); *B65D 25/205* (2013.01); *B65D*  
*45/16* (2013.01); *B65B 3/003* (2013.01); *B65B*  
*7/161* (2013.01)(21) Appl. No.: **15/428,045**(22) Filed: **Feb. 8, 2017****Related U.S. Application Data**(63) Continuation-in-part of application No. 15/425,863,  
filed on Feb. 6, 2017.(60) Provisional application No. 62/291,433, filed on Feb.  
4, 2016, provisional application No. 62/292,802, filed  
on Feb. 8, 2016.**Publication Classification**(51) **Int. Cl.***A61B 5/00* (2006.01)*B65B 7/16* (2006.01)(57) **ABSTRACT**

A method for producing prick test single use sterile vials for use in a prick test is provided. The method comprises providing at least one allergen source including therewithin allergen extract, providing a plurality of single dose vials, connecting a tube to the at least one allergen source, connecting a metering device to the tube, drawing by the tube an amount of the allergen extract from one of the at least one allergen source, receiving the amount of allergen extract at the metering device, and dispensing by the metering device a volume of allergen extract into one of the plurality of single dose vials.






FIG. 1




FIG. 2




FIG. 3




FIG. 4




FIG. 5




FIG. 6




FIG. 7




FIG. 8






FIG. 9








FIG. 11




FIG. 13




FIG. 14




FIG. 15




FIG. 16




FIG. 17




FIG. 18




FIG. 19




FIG. 20




FIG. 21A




FIG. 22






FIG. 21C




FIG. 21B




FIG. 23




FIG. 24







FIG. 25






FIG. 29




FIG. 30

## PRICK TEST SINGLE USE STERILE VIAL AND METHOD

### CROSS-REFERENCE TO RELATED APPLICATIONS

[0001] This application is a continuation-in-part of U.S. patent application Ser. No. 15/425,863, filed on Feb. 6, 2017, entitled ANTIGEN REGIONAL TESTING KIT. U.S. patent application Ser. No. 15/425,863 claims the benefit of and/or priority to U.S. Provisional Application No. 62/291,433, filed on Feb. 4, 2016, entitled ANTIGEN REGIONAL TESTING KIT. This application also claims the benefit of and/or priority to U.S. Provisional Application No. 62/292,802, filed on Feb. 8, 2016, entitled PRICK TEST SINGLE USE STERILE VIAL AND METHOD. Application Nos. 15/425,863, 62/291,433, and 62/292,802 are herein incorporated by reference in their entirety.

### TECHNICAL FIELD

[0002] The following disclosure relates to a sterile, single use prick test vial and method.

### BACKGROUND

[0003] Immunotherapy (IT) is recognized as one of the most curative treatment for allergies. By exposing the immune system to slowly increasing concentrations of immunomodulators such as an allergen or antigen, it will eventually stabilize and regain control the portion that is hypersensitive to the allergen or antigen. In general, immunotherapy is the “treatment of disease by inducing, enhancing, or suppressing an immune response.” Immunotherapies designed to elicit or amplify an immune response are classified as activation immunotherapies, while immunotherapies that reduce or suppress are classified as suppression immunotherapies. The active agents of immunotherapy are collectively called immunomodulators. They are a diverse array of recombinant, synthetic and natural preparations, often cytokines.

[0004] Immunotherapy involved in the treatment of allergies is a type of suppression immunotherapy, often termed desensitization or hypo-sensitization. This is compared with allergy treatments such as antihistamines or corticosteroids which treat only the symptoms of allergic disease. Immunotherapy is the only available treatment that can modify the natural course of the allergies, by reducing sensitivity to the immunomodulators such as antigens or allergens. An antigen and an allergen can both cause one's immune system to respond. An allergen is an antigen, but not all antigens are allergens. An antigen is any substance that is capable of causing one's immune system to produce antibodies. They are typically organic, or living, produced proteins. An allergen is any antigen that causes an allergic reaction. A non-allergen antigen could be a bacteria, virus, parasite, or fungus that causes an infection. This could also be something else that causes antibody immune system response, like toxins, chemicals, tissue cells involved in transplants or blood cells from a blood transfusion. An allergen is an environmentally produced substance that causes an allergic reaction, although the substance may not be harmful. Allergens cause no reactions in some individuals, while possibly causing a hypersensitive reaction in others. Common allergens include such things as pollen, plants, smoke, feathers,

perfumes, dust mites, toxic mold, food, drugs, animal dander, and insect bites and stings.

[0005] The exact mechanisms of how IT works are not fully understood, but they involve shifting a patient's immune response from a predominantly “allergic” T-lymphocyte response to a “non-allergic” T-lymphocyte response.

[0006] Current accepted processes for performing allergy immunotherapy include injecting immunomodulator matter in the form of antigen material into patient subjects. This is referred to as subcutaneous immunotherapy (SCIT), requiring a patient to visit a doctor's office for weekly injections. It is very expensive and time-consuming. A second technique, sublingual immunotherapy (SLIT), involves the application of allergy extracts (antigens), and allergens placed into a pill form and swallowed by the patient or disposed in “allergy drops” which are placed under the tongue for the allergens/antigens to be absorbed into the oral mucosa. Transdermal patches may have been used without much success and mostly were used for patch testing to see if a patient reacts to various chemicals or allergens.

[0007] Of the people who start traditional subcutaneous injected immunotherapy (SCIT), 90% fail to complete their therapy due to needle fatigue and not being able to see a doctor in their office once or more per week for several years. Further, doctors charge for every one of those visits. Further, doctors trained to give injections for allergies are concentrated in high population and upper middle class places. People in rural areas and people who do not live in upper middle class areas cannot get to an allergist for shots. Consider an inner city kid having to ride public transportation and pay a high copay just to get a high risk injection if an alternative therapy were available!

[0008] Allergies are also linked to depression and suicide and are among the top ten reasons for missed work and lost productivity. Lastly, allergies and asthma result in billions of dollars in lost productivity and healthcare costs among the 90% of allergy patients that either never get immunotherapy or fail immunotherapy delivered under its current administration methods.

[0009] Skin allergy testing is a method for medical diagnosis of allergies that attempts to provoke a small, controlled, allergic response. A microscopic amount of allergen is introduced to a patient's skin by various means. One of these is a prick test wherein an allergen is introduced by pricking the skin with a needle to introduce a small amount of allergen into the epithelial layer of the skin. Similar tests include the scratch test that involves a deep dermic scratch made with the blunt bottom of a lancet, an intra-dermic test involving a tiny quantity of allergen injected under the dermis with a hypodermic syringe, a scrap test involving a superficial scrape performed with the help of the bevel of a needle to remove the superficial layer of the epidermis and a patch test wherein a patch is applied to the skin, the patch containing a small amount of allergen.

[0010] In the prick test, a small amount of purified allergen is generally pricked onto the skin surface, usually the forearm. This test is usually done in order to identify allergies caused by such things as pet dander, dust, pollen, foods or dust mites. The skin prick test involves first placing a small amount of substances that may be causing the symptoms on the skin, most often on the forearm or upper arm back. Some of the prick test kits that are available are utilized to maximize the number of tests that can be per-

formed with a given amount of allergen. Thus, a 5 mL bottle of allergen can provide up to 1800 tests. Typically, a number of wells formed in a tray are each filled with an allergen, for example with about 0.125 mL of extract. The tip of a pick or lancet is dipped into the small amount of extract disposed in the bottom of the well and then a particular site is pricked on an individual. In some cases multiple pricks are mounted in or on a frame to provide for multiple well associations, such that a series of six or seven wells in a line in the tray can be simultaneously accessed so that a number of different allergens can be applied to an individual's skin simultaneously. There are even larger arrays of picks that can be dipped into the wells for treating a large area with more different allergens at the same time. After use, the wells are typically sealed by inserting a prick into the well.

**[0011]** One of the issues involved in prick testing is that the allergen must be removed from a sterile bottle with a syringe to fill each well. The sterile bottle containing the allergen typically has some type of injection stopper, typically rubber, through which the hypodermic needle is inserted to extract the allergen. Once the hypodermic needle is inserted through the injection stopper, the sterile barrier is broken. The allergen is disposed into the well and the pick disposed therein multiple times. There is not necessarily any sterilization procedure between applying the allergen to the skin of an individual and then disposing the pick back into the well to seal that well. Thus, once the allergen is extracted from the 5 mL bottle, the first use of the allergen is the only sterile use. Thereafter, application thereof of the allergen up to **1800** times can potentially result in issues.

#### SUMMARY

**[0012]** In one aspect thereof, a method for producing prick test single use sterile vials for use in a prick test is provided. The method comprises providing at least one allergen source including therewithin allergen extract, providing a plurality of single dose vials, connecting a tube to the at least one allergen source, connecting a metering device to the tube, drawing by the tube an amount of the allergen extract from one of the at least one allergen source, receiving the amount of allergen extract at the metering device, and dispensing by the metering device a volume of allergen extract into one of the plurality of single dose vials.

**[0013]** In another embodiment, the method further comprises sealing the one of the plurality of single dose vials with one of a plurality of stoppers and securing the one of the plurality of stoppers on the one of the plurality of single dose vials with one of a plurality of seals.

**[0014]** In another embodiment, the method further comprises affixing one of a plurality of labels to the one of the plurality of single dose vials.

**[0015]** In another embodiment, the method further comprises providing a tray having disposed therein a plurality of recesses and placing within one of the plurality of recesses of the tray the one of the plurality of single dose vials.

**[0016]** In another embodiment, the method further comprises repeating the drawing, receiving, dispensing, sealing, securing, affixing, and placing steps until a desired number of the plurality of single dose vials are produced.

**[0017]** In another embodiment, the method further comprises providing a prick tester having a plurality of needles extending thereon, aligning the plurality of needles of the prick tester with the plurality of single dose vials, inserting each of the plurality of needles of the prick tester into the

plurality of single dose vials, and applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

**[0018]** In another embodiment, the method further comprises extracting a predetermined amount of allergen extract from one of the plurality of single dose vials and disposing the predetermined amount of allergen extract into a corresponding one of a plurality of wells.

**[0019]** In another embodiment, the method further comprises repeating the extracting and disposing steps until a desired number of the plurality of wells contain allergen extract.

**[0020]** In another embodiment, the method further comprises providing a prick tester having a plurality of needles extending thereon, aligning the plurality of needles of the prick tester with the plurality of wells, inserting each of the plurality of needles of the prick tester into one of the plurality of wells, and applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

**[0021]** In another embodiment, the method further comprises providing a prick tester having a plurality of needles extending thereon, aligning the plurality of needles of the prick tester with the plurality of single dose vials, inserting each of the plurality of needles of the prick tester into the plurality of single dose vials, and applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

**[0022]** In another embodiment, the method further comprises extracting a predetermined amount of allergen extract from one of the plurality of single dose vials and disposing the predetermined amount of allergen extract into a corresponding one of a plurality of wells.

**[0023]** In another embodiment, the method further comprises repeating the extracting and disposing steps until a desired number of the plurality of wells contain allergen extract.

**[0024]** In another embodiment, the method further comprises providing a prick tester having a plurality of needles extending thereon, aligning the plurality of needles of the prick tester with the plurality of wells, inserting each of the plurality of needles of the prick tester into one of the plurality of wells, and applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

**[0025]** In another aspect thereof, a system for producing prick test single use sterile vials for use in a prick test is provided. The system comprises an allergen source including therewithin allergen extract, a tube connected to the allergen source to draw the allergen extract from the allergen source, a metering device connected to the tube, wherein the metering device receives by way of the tube an amount of allergen extract from the allergen extract in allergen source, and wherein the metering device serves to dispense a volume of extract, a plurality of single dose vials, wherein each one of the plurality of single dose vials is configured to receive the volume of extract dispensed by the metering device, a plurality of sterile injection stoppers for use in sealing the plurality of single dose vials, a plurality of seals to secure the plurality of sterile injection stoppers on the plurality of single dose vials, a plurality of labels, each one of the plurality of labels configured to be affixed to one of the plurality of single dose vials, a tray having disposed therein a plurality of recesses, each one of the plurality of

recesses able to receive one of the plurality of single dose vials, and a prick tester having a plurality of needles extending thereon, wherein the plurality of needles are arranged to be inserted into the plurality of single dose vials disposed in the plurality of recesses of the tray.

**[0026]** In another embodiment, the system further comprises a dilutant for dispensing into a number of the plurality of single dose vials.

**[0027]** In another embodiment, the dilutant comprises a saline solution.

**[0028]** In another embodiment, the dilutant comprises a glycerol solution.

**[0029]** In another embodiment, the system further comprises a first control material and a second control material, wherein the first control material is dispensed into a number of the plurality of single dose vials, and wherein the second control material is dispensed into a number of the plurality of single dose vials.

**[0030]** In another embodiment, the first control material is water.

**[0031]** In another embodiment, the second control material is a histamine.

#### BRIEF DESCRIPTION OF THE DRAWINGS

**[0032]** For a more complete understanding, reference is now made to the following description taken in conjunction with the accompanying Drawings in which:

**[0033]** For a more complete understanding, reference is now made to the following description taken in conjunction with the accompanying Drawings in which:

**[0034]** FIG. 1 illustrates a diagrammatic view of a dilution sequence of diluting a concentrated antigen extract;

**[0035]** FIG. 2 illustrates a diagrammatic view of a production line for filling distribution bottles;

**[0036]** FIG. 3 illustrates a flow chart for the dilution process;

**[0037]** FIG. 4 illustrates a flow chart for the dispensing process;

**[0038]** FIG. 5 illustrates a diagrammatic view of a color-coded box with different diluted levels of allergens;

**[0039]** FIG. 6 illustrates a cross-sectional view of a low dose bottle;

**[0040]** FIG. 7 illustrates a diagrammatic view of filling sterilized bottles with concentrated antigen;

**[0041]** FIG. 8 illustrates a top view of a testing kit illustrating the concentrated antigen and the prick test wells;

**[0042]** FIG. 9 illustrates a diagrammatic view of the operation of transferring concentrated antigen from one of the concentrated antigen bottles in the kit to an associated well in the kit;

**[0043]** FIG. 10 illustrates a side cross-sectional view of the wells after antigen is disposed therein and illustrating applying the antigen in the wells to the tips of a prick tester;

**[0044]** FIG. 11 illustrates a perspective view of the wells after filling in the disposal of a prick tester over a select ones of the wells;

**[0045]** FIG. 12. Illustrates a top view of the embodiment of FIG. 11;

**[0046]** FIG. 13 illustrates a top view of a well portion of the kit with the protesters disposed thereon for storage;

**[0047]** FIG. 14 illustrates a side view of the overall well portion and the stored prick tester with a sterile seal disposed thereon;

**[0048]** FIG. 15 illustrates a diagrammatic flow chart of the overall operation of assembling the kit and utilizing the kit;

**[0049]** FIG. 16 illustrates a diagrammatic view of the kit and a close configuration;

**[0050]** FIG. 17 illustrates a detailed view of a drawer for the information sheet;

**[0051]** FIG. 18 illustrates an exploded view of the kit with the addition of a dilutant addendum;

**[0052]** FIG. 19 illustrates a side view of the fully closed kit;

**[0053]** FIG. 20 illustrates the dilution process in a step-wise sequence;

**[0054]** FIGS. 21A-21C illustrate top views of different types of kits for standardized antigens and region specific antigens and the various assemblies;

**[0055]** FIG. 22 illustrates a detailed view of the various regional kits;

**[0056]** FIG. 23 illustrates a mapping view of the United States with the various allergy regions;

**[0057]** FIG. 24 is a diagrammatic view illustrating the fill operation for producing a sterile, single use vial of allergen for prick testing;

**[0058]** FIG. 25 is a diagrammatic view of a first seal operation for producing a sterile, single use vial of allergen for prick testing;

**[0059]** FIG. 26 is a diagrammatic view of a stoppering operation for producing a sterile, single use vial of allergen for prick testing;

**[0060]** FIG. 27 is a diagrammatic view of a labeling operation for producing a sterile, single use vial of allergen for prick testing;

**[0061]** FIG. 28 is a diagrammatic view of a packaging operation for producing a sterile, single use vials of allergen for prick testing;

**[0062]** FIG. 29 is a flow diagram illustrating one method of producing a sterile, single use vial of allergen for prick testing; and

**[0063]** FIG. 30 illustrates a filled, stoppered, and labeled sterile, single use vial of allergen for prick testing.

#### DETAILED DESCRIPTION

**[0064]** Referring now to the drawings, wherein like reference numbers are used herein to designate like elements throughout, the various views and embodiments of an antigen regional testing kit are illustrated and described, and other possible embodiments are described. The figures are not necessarily drawn to scale, and in some instances the drawings have been exaggerated and/or simplified in places for illustrative purposes only. One of ordinary skill in the art will appreciate the many possible applications and variations based on the following examples of possible embodiments.

**[0065]** Referring now to FIG. 1, there is illustrated a depiction of a typical technique for diluting immunomodulators such as antigens, as one example. Preparation of a diluted antigen is performed first by receiving a bottle of extract concentrates from an approved vendor. These are formulated in a given weight/volume (w/v) format with a given antigen associated therewith. For typical antigens such as those associated with the cat antigen, these are relatively well controlled. Typically, a vendor will provide an extract for a single antigen or allergen. Allergens such as pollen and the such are not as well controlled due to the technique for collecting such. In any event, there are typically very few approved vendors for these extracts and allergist typically



receives these vendor provided concentrates in a sufficient quantity to make the necessary diluted solution.

**[0066]** Allergen extract is typically comprised of a non-allergenic material, a non-allergenic protein, and an allergenic protein. The extraction solutions can be aqueous containing saline and phenol which could be a glycerinated solution. The allergen is added, the units of measure are sometimes referred to as "AU" for "allergy units," typically used for mites. These are referred to as "AU/mL." For such things as grass and cats, the term "BAU" is used for "bioequivalent units." For other allergens, the terminology is, for example, 1:20 w/v, which stands for 1 g source material per 20 mL of fluid. The relationship between BAU and 1:20 w/v depends upon the extract. In any event, there is a defined amount of extract contained within the concentrate.

**[0067]** When concentrated extracts are formulated by an authorized vendor, they are typically provided in standardized versions and non-standardized versions. In standardized versions, they typically are provided in a 50% glycerin dilutant. They can either be a single allergen extract, or they can be a mix. For example, one can obtain a "9 Southern Grass Mix (concentrate)" which contains equal parts of: 2 Bermuda at 10,000 BAU/mL, P 27 7 Grass at 100,000 BAU/mL, and 15 Johnson at 1:20 w/v. For non-standardized extracts, these are typically provided in either a glycerin dilutant or an aqueous dilutant such as saline. They can be a single extract or a mix. Thus, whenever a concentrated extract is referred to hereinbelow, this refers to a formulation that is provided by an authorized vendor that can be diluted in accordance with the processes described hereinbelow. These are typically provided in the 50 mL bottles with a needle compatible.

**[0068]** Referring back to FIG. 1, the extract concentrate is disposed in a bottle 102. This is a sterile concentrate that has an injection stoppered top 104. There are provided a plurality of five 5 mL sterile injection stoppered bottles 106, 108, 110, 112 and 114, although there could be more, and the bottles or containers could be larger than 5 mL. Each of these bottles has disposed therein a defined amount of dilutant, depending upon what the final dilutant is required to be. Typically, the amount of dilutant is 4.5 mL. The procedure is to, first, extract a defined amount of the concentrated extract from the bottle 102 and dispose it in the bottle 106. This is facilitated by the sterile hypodermic that is inserted through the stopper at the top of the bottle 102 to extract concentrate, and then the hypodermic is inserted to the stopper in the bottle 106 to inject extract from bottle 102 into bottle 106. Typically, the concentration in the concentrated extract bottle 102 is 1:20 w/v. This will result in a dilution of 1:10 in bottle 106. Then, 0.45 mL of the diluted solution from bottle 106 is extracted and inserted into bottle 108, resulting in a 1:100 dilution of the original concentrate in model 108. The process is repeated up to the bottle 114 to provide a solution that is at a dilution of 1:100,000 of the original concentrate. This is a conventional way to provide a selected dilution of the original antigen. However, it should be understood that any concentration level can be provided from one bottle to the next. The purpose of using the sequential bottles is to allow an achievable portion of one bottle to be distributed to the next bottle, rather than trying to extract a very small amount of the initial concentrated extract. Typically, an allergist will then extract from the desired dilution an amount of the diluted antigen for injection

percutaneously. Typically, desensitization is achieved by using the most diluted antigen level initially and sequentially moving up to a higher concentration level over time.

**[0069]** Illustrated in FIG. 1 are three hypodermic needles, one selecting a "dose" from bottle 114, and labeled hypodermic 116, a second hypodermic needle 118 for retrieving a dose from bottle 112, a third hypodermic needle 120 for extracting a dose from bottle 110. Each of the hypodermic needles 116, 118 and 120 will contain a different diluted dose. These would typically be separate needles in the event that the allergist or medical professional is injecting a patient. For other purposes, they could be the same needle, depending upon the dose or concentration required. A "dose" is defined by the amount of all the diluted product that would be required for the desired immunotherapy. This is defined by the medical professional. If, for example, bottle 112 were utilized, it may be that 1 mL of diluted solution constituted a "dose." It could be that less than 1 mL constituted a "dose."

**[0070]** In general, the typical distribution chain requires that the allergist or other medical professional purchase the base concentrate and then perform the dilution process. However, this procedure typically requires breaking the seal on the base concentrate bottle and then inserting a needle into the base concentrate bottle for the first dilution step. This occurs multiple times. Thus, multiple needles, each being sterile, can be used one time or, more commonly, a single needle is utilized in association with the base concentrate bottle, with the assumption that, since it does not involve insertion into human flesh, it is still sterile. In any event, this needle must penetrate the rubber stopper seal on the base concentrate bottles multiple times. In fact, these bottles could typically be held upside down and they would leak and, once the seal is broken, there is no sterile cover over the rubber stopper. This is a result of the multiple needle piercings of the rubber stopper. This is also the case with the small 5 mL bottles in that each has to be penetrated at least twice in the higher concentrate bottles. Thus, the last bottle that the allergist has would be a 5 mL bottle, and this bottle would already have one piercing of the rubber stopper seal in order to provide the initial dilution level into the carrier material, such as saline. Thereafter, a patient might be able to receive 5 or 10 doses from that particular bottle, requiring 5 or 10 more piercings of the rubber stopper. During this time, of course, there is no seal over the rubber stopper.

**[0071]** In order to solve this problem, a process is provided whereby the concentrate bottles are produced from a sterile environment container which already has the dilutant provided thereto. This is illustrated in FIG. 2. In this embodiment, a large bottle of base concentrate 202 is provided at a first concentration. This is then metered into a first dilutant reservoir 204 via a tube 206 and a metering valve 208. The first dilutant container 204 has contained therein a carrier liquid 210 which can be, as described herein below, any of a number of materials, such as saline, glycol and the such. Typically, this depends upon the base carrier material associated with the concentrate 202. By knowing the volume of the material to 10 within the reservoir 204 prior to metering therein of the concentrate 202 with the metering valve 208, a very accurate amount of concentrate 202 can be dispensed within the reservoir 204. This metering valve 208 can be controlled to such a level that a very fine and controlled concentrated level can be defined. At this point in the

process, the concentrate **202** defines a “batch” material that, for some allergens, is important. For example, if the allergen was related to pollen or the such, this can vary depending upon the year, the production harvest, the quality of harvest, etc. By defining a batch, and controlling the quality and the concentration level at each step in the dilution level, a very controlled dilution level can be provided for that particular batch.

**[0072]** Once the concentrated level or dilution level in the bottle has been defined, this is then utilized to provide a controlled amount of diluted allergen to a second dilution bottle **216** through a metering valve **218**. This is repeated for multiple bottles down the line to a last bottle **220**. Thus, there are then provided a plurality of larger vessels with controlled dilution levels at each diluted stage in a sterile environment. Each of these bottles **204**, **216**, and **220** has associated therewith a control metering valve **222** and a dispensing nozzle **224** that is operable to dispense diluted allergen material into a receptacle. This metering valve **222** and associated nozzle **224** are all approved to interface with an approved bottle.

**[0073]** In the dispensing process, there are provided for each dilution stage a plurality of bottles **230**. Each of these bottles has a shape and opening that is approved to be interfaced with the nozzle. Each is passed by the nozzle and an exact amount of diluted material, the allergen, dispense therein. Thereafter, each bottle **230** is then subjected to a Procedure which inserts the rubber stopper and the seal there over. The result is a bottle with a defined dosage amount at a defined concentration level for a known batch of allergen, all of which is a sterile environment. The allergist or medical professional need only then remove this seal, insert the needle through the rubber stopper and extract the appropriate amount of diluted allergen. This particular bottle **230** is designed to be a single dose. Thus, the actual amount of material disposed within the bottle will be approximately 1 mL. Typically, the bottles are 5 mL bottles. In this situation, that 5 mL bottle can be modified to maintain the same size but only provide for a single 1 mL dose. This will be described hereinbelow.

**[0074]** Referring now to FIG. 3, there is illustrated a flowchart for the overall dilution process, which is initiated at a block **302**. The process flows to a block **304** to provide the base concentrate. As noted hereinabove, the base concentrate is typically associated with a fixed batch so that it is a well-controlled product, and this is a sterile bottle within a sterile environment, and the size can be larger than typically distributed to an allergist. The process flows to a block **306** wherein the various dilutant containers are provided with a sterilized base liquid such a saline or glycol, depending upon intended final use. The process then flows to block **308** wherein a metered amount of dilutant is passed from the base concentrate vessel to the first dilutant container. The process then flows to a block **310** wherein a metered amount of diluted allergen is transferred from the previous bottle to the next bottle until the last bottle is complete, which is determined at a decision block **312**. This will then result in multiple vessels of diluted allergen at the appropriate diluted level for final dispensing. The program then flows to a block **316**.

**[0075]** Referring now to FIG. 4, there is illustrated a flowchart for the dispensing operation, which is initiated at a block **402**. The process flows to a block **404** wherein a nozzle is inserted from the select dilutant container or vessel

to the final dispensing bottle, typically a single dose bottle. The process then flows to a block **406** to meter a finite dose to the end user bottle. Again, this could be a 1 mL dose in a 5 mL bottle or a 1 mL dose in a 1 mL bottle. The process then flows to a block **408** to send the end user bottle to capping, followed by a process block **410** to dispose a cap on the bottle, this being a rubber seal to allow insertion of the needle there through, and then to a seal block **412** to provide a sterile seal over the bottle and then to a terminate block **416**.

**[0076]** Referring now to FIG. 5, there is illustrated a depiction of a delivery container for the final allergen in the bottles **230**. The bottles would typically be controlled such that they would have, for example, a very distinctive color associated with each diluted level. The highest concentration bottle would be a bottle that would have to be carefully dispensed to a patient, as, if not adequately desensitized to the allergen, this could result in anaphylactic shock to the patient. Thus, it is important that the correct dose at the correct dilution level is administered at a particular time. This color coding process for any type of identification process clearly presents to the medical professional some clear indication of the concentrated level of allergen. Typically, these allergens will be provided in a kit form. For example, it may be that a patient would require the lowest concentrate level of the allergen three times per week for two weeks, followed by the next concentrated level of allergen two times per week for one week, etc. Thus, the various diluted levels of concentrate would be provided in single-dose bottles to, for example, a pharmacy, which would dispense the particular dosages. The first thing is that they would be from a common batch and they would be provided to the patient in the appropriate presentation. For example, the patient might receive a first box of six bottles **230** at the lowest concentrate level for the first two weeks. The pharmacist then would provide the second part of the prescription for the second level in the form of two bottles of concentrate at the next level, this all being color-coded. The box might be color-coded, as well as the bottles. Again, each of these bottles would be a single dose. Since they are all single-dose and contained within sterile bottles, the shelf life is considerably longer.

**[0077]** Referring now to FIG. 6, there is illustrated a cross-section of a typical bottle that would be involved with respect to providing a single dose of 1 mL in a larger standard 5 mL bottle. A 5 mL bottle is represented by a bottle **602**. This bottle **602** has an opening **604** that is provided for the 5 mL bottle. An insert **606** is provided for filling space within the bottle. This can be any type of plastic insert, etc. A smaller insert bottle **608** is disposed within the insert **606** to provide an elongated interior **620** which has a volume slightly in excess of 1 mL, such that a 1 mL dose can be disposed therein. This elongated interior is covered with a rubber stopper **622** such that a needle **624** can be disposed there through and be able to extract 1 mL of diluted antigen. If not for the elongated opening **620** facilitated by the insert **606** in the bottle **608**, the 1 mL of diluted antigen would be disposed at a lower level and would be more difficult to extract.

**[0078]** FIG. 7 illustrates a diagrammatic view for filling sterilized bottles with concentrated antigen. For a particular antigen, a vessel **702** is provided containing the concentrated antigen which is fed through a tube **704** via a cut off **706** and a cutoff **708** to a dispensing nozzle **710**, similar to the

dispensing nozzles **224** described hereinabove. A plurality of empty sterilized bottles **712** that are passed along an assembly line underneath the nozzle **710**, which dispenses a predetermined amount of concentrated antigen therein. This results in a plurality of the bottles **712** with a dispensed amount **714** of antigen disposed therein. These bottles **712** can range in size from a 0.5 mL bottle to a 10 mL bottle. The preferred size is approximately 5 mL. This, of course, depends upon the application, as will be described hereinbelow.

[0079] The next step in the process is to bottle **712** with a rubberized stopper **716**, as described hereinabove with respect to FIG. 2. This results in a sterilized bottle **720**, noting that the entire process is carried out in a sterile environment. At this point, the rubberized stopper **716** on the filled bottle **720** is covered with a layer foil to basically provide an overall sterile vessel. As noted hereinabove, to use this bottle **720** after sterilization, all that is required is to peel away the foil top and insert the needle through the rubberized stopper **716**. At this point, of course, the bottle **720** will be in an unsterile environment. Thus, the first puncture with a sterile needle into the rubberized stopper **716** will be a sterile operation. Thereafter, a sterile environment cannot be ensured. At this point, as will be described hereinbelow, the bottle **720** is labeled or color-coded. Since this bottle **720** will be disposed in a kit, is important that there be some indicia provided that will allow the antigen contents of the bottle to be discriminated from the contents of the other bottles with different antigens disposed therein. This indicia can be an alpha character, a numeric character, a color reader barcode.

[0080] Referring now to FIG. 8, there is illustrated a top view of a clamshell configuration having two portions, a portion **802** and a portion **804**. These are joined together by a hinge **806** such that the clamshell configuration can be closed. The portion **802** is operable to contain a plurality of discrete antigens in bottles **720**. There are also provided two bottles **810** that are control bottles. One of the bottles will contain water and the other will contain a histamine. The water, as the control, applied to the skin should elicit no response whereas the histamine, when applied to the skin, will elicit a guaranteed response. This provides the medical technologist some type of benchmark for the particular individual.

[0081] The other side **804** constitutes a portion that contains a plurality of wells **812** that have disposed therein a very small portion, a drop or two, of antigen, wherein each of the wells **812** is associated with the corresponding bottle **720** with an antigen. Above each well **812** is provided an indicia marker **814** which is associated with the indicia on each of the bottles. It can be that the indicia is provided not only on the bottle but also on the insert associated with the bottle **720** in the portion **802**. Thus, as illustrated in FIG. 9, a particular bottle **720** is retrieved from a well **820**. Again, this well could have an indicia **822** disposed thereabove that is associated with a particular antigen for a particular kit. A hypodermic **824** is utilized to extract a portion of the antigen disposed within the bottle **720** through the rubberized stopper **716** after the foil is removed therefrom, and this extracted antigen is disposed within a corresponding well **812**. Again, with color coding, it is possible to have the bottle coded for a certain color, the receiving pocket **820** color-coded for that bottle, and the well **812** associated therewith color-coded for that antigen. In such a manner,

when the antigens are assembled for a particular kit, it is possible to ensure that the appropriate antigen is inserted in the appropriate pocket **820** by matching the colors, and, when utilizing the antigen to fill the well for a prick test, it can be ensured that the appropriate color coding is matched, since the bottle itself will have color coding disposed thereon. Again, an alpha character or numeric character could be utilized in the same manner.

[0082] As can be seen from FIG. 9, all that is required for the medical technologist is to open up a particular clamshell prick test kit and utilize the appropriate concentrated antigen for dispensing in the well **812**, this dispensing requiring 1, 2 or 3 drops—a very small amount. Thereafter, as illustrated in FIG. 10, a prick tester **1002** is provided having a plurality of needles **1004** extending thereon and aligned with at least a row of wells **812**. This protector **1002** is disposed there over and the edges aligned with alignment slots **1005** and **1006**, there being a corresponding alignment tab **1008** and **1010** on the prick tester **1002** to ensure that needles **1004** are correctly aligned with the center of the well **812**. When the prick tester **1002** is disposed down over the portion **804** and the tabs **1008** and **1010** are aligned and inserted within the slots **1005** and **1006**, respectively, the tips of the needles **1004** will be sufficiently disposed within the corresponding well **812** in order to ensure that the tip of the needle has concentrated antigen disposed thereon. This is then applied to the skin of the patient to elicit the potential response thereto.

[0083] Referring now to FIG. 11, there is illustrated a perspective view of the portion **804** having associated therewith the wells. In this portion, each of the wells **812** are disposed in rows, there being four rows each having four wells **812**. Below those four rows are provided a single row for the two control wells **813**. There is illustrated one configuration of the prick tester which allows a single row of wells **812** to be accessed. In this embodiment, a prick tester **1102** is provided that only accesses one row of wells and is disposed thereover with, in an exemplary embodiment, an indicia of an alpha character, A, B, C and D. Disposed thereunder will be a corresponding well **812** with the corresponding alpha character. The tab **1005** will align with the corresponding slot in the portion **804** associated with that row. It is noted that there will be provided a separate prick tester **1102** for each of the rows with corresponding matching alpha characters. In order to ensure that the prick tester will only be insertable into the portion **804** for the associated row, different dimensioned slots **1104** will be provided on one side of the row and different dimensioned slots **1106** will be disposed on the opposite side of the row. Thus, the prick tester **1102** for each row will not only be designed to only be associated with that row, but it also has to be oriented in the correct direction such that the rows are correctly designated for a particular prick tester. It could be that the prick test **1102** could cover two rows or could cover all four rows or even all four rows **812** and the control wells **813** in the bottom row.

[0084] Referring now to FIG. 12, there is illustrated a top view of the portion **804** illustrating not only the wells **812** but also a storage region **1202** for each of the rows beneath the particular row. It will be a corresponding slot disposed in this particular storage region **1202**, that storage region **1202** being a row in and of itself. Thus, for example, there will be a slot **1005** disposed on one side of the top row and a slot **1006** disposed on the opposite side of each of the wells

**812** in that row. Directly therebeneath, there will be a corresponding slot **1005'** on one side of the storage region **1202** associated with that first row, and a corresponding slot **1006'** disposed on the opposite side. There will be disposed directly beneath each well **812** a receiving hole **1208** for receiving the needle on the bottom of the associated prick tester.

[0085] Referring now to FIG. 13, there is illustrated a top view of the portion **804** in a shipping configuration. In this configuration, there will be provided a prick tester **1302** for the first row, a prick tester **1304** for the second row, a prick tester **1306** for the third row, a prick tester **1308** for the fourth row and a prick tester **1310** for the control row. Each of these is disposed in the respective storage region **1202** associated with each particular row. It is only necessary thereafter to remove the prick tester **1302** after concentrated antigen has been disposed within all or select ones of the wells in the first row to allow utilization thereof.

[0086] Referring now to FIG. 14, there is illustrated a side view of the portion **804** in the shipping configuration with the prick tester **1310** illustrated from the side. This is assembled in a sterile environment and, thereafter, it is sealed with a see-through seal **1402** that maintains the sterile environment. This can be removed when it is being used. Of course, there is no seal required over the antigen bottles **720** in the portion **802**, as they are already sealed individually.

[0087] Referring now to FIG. 15, there is illustrated a process flow illustrating the process of dispensing the concentrated antigen into the bottles **712** and then sealing them to provide the final bottle **720** with the rubberized stopper **716** disposed thereon. These are then each marked with the appropriate indicia. There can be a procedure wherein each antigen has a unique color associated there with some type of unique code. Due to the large number of antigens, it may be that two colors would have to be provided to provide discrimination among a large number of antigens. However, for standardized kits that require a finite number of antigens, this may not be necessary. If a large number of kits of a particular regional nature are assembled, it can be that the concentrated antigens are color-coded for those kits and maintained in an identified area. They can then be assembled into the kits as necessary. This is illustrated with directional arrows indicating that, for example, two wells **1502** and **1504** have bottles **720** disposed therein associated with antigen A and antigen B, respectively, for example. These are not necessarily directly disposed therein after filling but, rather, could be filled from a queuing station. The control bottles would, of course, be fairly standardized and may even have a standardized color across all kits. The next step is to then extract antigen from the select one of the bottles **720** and dispose it within a corresponding well **812**. The prick tester **1002** and the pin **1004** associated with that well **812** will then be mated to coat the tip of needle **1004** with the small amount of antigen that was disposed in the well **812**. This is then applied to the skin **1510** of a patient.

[0088] Referring now to FIG. 16, there is illustrated a closed view of the clamshell configuration of the two portions **802** and **804**. They will be disposed such that they will come together and be slightly separated with tabs **1602**. A securing latch **1604** is disposed alongside to secure the two portions **802** and **804** together facing each other. There will be a label **1610** disposed on the upper side thereof identifying the particular point. In addition, there will be provided, in one example, and information slot **1612** within

the portion **802**, although it could be in the portion **804** or disposed any place with respect to the overall clamshell configuration. A detail of this information slot **1612** is illustrated in FIG. 17. This is illustrated as a drawer **1702** that can be pulled out with an information brochure **1704** disposed therein. This information brochure could contain advertisements, instructions and also well maps for the various wells indicating the antigens that are disposed therein for this particular kit. There are also be provided safety instructions and emergency numbers and the such.

[0089] Referring now to FIG. 18, there is illustrated a view of an open clamshell with an addendum section **1802** attached thereto. This addendum section **1802** has disposed therein a dilutant set of saline bottles **1806**, although they could contain other dilutants that are useful for the dilution operation, as described hereinabove. This is hinged to the side of the portion **804** by a hinge **1810**. Additionally, the information section is disposed in a section attached by a hinge **1814** to a section **816** that is foldable. When this entire system is folded into a multiple level clamshell configuration, the configuration is as illustrated in FIG. 19.

[0090] FIG. 20 illustrates a diagrammatic view of the operation of creating a set of diluted antigens in accordance with a therapeutic regimen provided as a prescription. In accordance with the operation of the kit, the first step is to determine the sensitivity of a particular individual to one or more of the antigens disposed in the kit. This is the "testing" phase. Once it is determined which antigens elicit a response, a cocktail can then be created. This cocktail of antigens is facilitated by selecting the first bottle **1806** from section **1802**. As noted in FIG. 18, each of the bottles **1806** is disposed within a well **1820**. Above each of these wells is disposed a designator D1-D6. These are associated with the dilution level of a particular model. Each of these bottles has contained therein a predetermined amount of saline solution. The typical procedure, as described hereinabove, is to first mix concentrated antigen within the first bottle **1806** from select ones or more of the antigen bottles **720**. This is facilitated with a hypodermic needle. The same hypodermic needle could be utilized for each extraction, but there may be some issue with sterility and a separate hypodermic needle could be utilized for each extraction. This is up to the medical technologist. A predetermined amount of concentrated antigen is extracted from the associated antigen bottle **720** that elicited a response during the prick test in the testing phase and is then inserted into the bottle **1806** through the associated stopper **1807** on each bottle **1806**. The bottle D1 is initially filled with approximately 8 mL of saline solution or similar dilutant. The bottle is a 10 mL bottle. To achieve the appropriate dilution level in the bottle **1806**, the instructions will inform the medical technologist how much dilutant to extract from each of the bottles. In one example, a level of 0.2 mL could be extracted, allowing 10 extractions to be inserted within the bottle **1806**. If all 12 bottles were utilized, this could result in 2.4 mL of concentrated antigen from 12 bottles **720** be inserted into the bottle **1806**. Thus the bottle **1806** could be filter level of 7.6 mL of dilutant initially. This resulted in an initial level of 7.6 mL filled potentially to 10 mL. A finite amount of dilutant is extracted from the bottle **1806** labeled D1 and transferred to the second bottle **1806** labeled D2 and so on until the last transfer to the bottle **1806** labeled D6. This will be a final bottle. These bottles could also be color-coded or marked with the particular dilutant level such that they could be

placed back into the kit. The reason for this is that a therapeutic program would require a patient to initially receive injections from the highest diluted bottle **1806**, associated with the level D6 in three or four doses or more of 0.5 mL each. The patient would then move on to the higher next higher diluted level and so on. This was described hereinabove. Thus, to recap, the procedure is to provide a test kit to a pharmacy which could then be distributed to a medical technologist. This kit can provide all of the antigens necessary to perform the prick test and the tools to facilitate that testing operation in the testing phase. If this was just a prick test kit and did not require any dilution process, only a small bottle would be required, on the order of possibly 0.5 mL or 1 mL. Thus, the bottle **720** could be very small. However, by providing the dilutant section **1802** in the clamshell with the dilutant bottles, larger bottle **720** could be facilitated to include a higher level of antigen for the purpose of creating a dilution set. Again, these could be as high as 5 mL or 10 mL. In some situations, a 10 mL bottle **720** may be facilitated really due to the fact that that is the industry standard. However, the size of the bottle merely reflects the amount required to, first, perform the prick test and, second, create the dilution sets. Since only 0.2 mL is required to create the dilution sets from any one bottle **720**, it may be that a small bottle of less than 1 mL could be utilized. The clamshell configuration therefore provides an entire solution to either just the prick test or a combination of the prick test and the dilution set creation. Thus, a diagnosis can be performed and then a prescription created and filled all with a single kit.

**[0091]** Referring now to FIG. 21A, there is illustrated a top view of a first section **802** having disposed to sections therein. There are provided eight wells **820** associated with a standardized section **2102** and 7 wells **820** associated with a region specific section **2104**. The reason for this is that the kits can be “regionalized.” This will allow a particular kit to be provided to a given region, wherein standardized antigens could be provided which are utilized over all regions and region specific antigens can be provided for a particular region and even a particular season in that region. The controls will be common across all kits. FIG. 21B illustrates a second configuration which is just standardized antigens. It can be seen less than 12 antigens might be included in this kit, as illustrated in un-dashed wells **820**. Similarly, there can be provided a regional kit illustrated in FIG. 21C wherein only the wells **820** that are un-dashed are filled.

**[0092]** When the kits are originally assembled, each of antigens has an NDC associated therewith. When all of the antigens are assembled into a “testing” kit, a separate NDC is associated with that kit. This particular NDC associated with that kit is an approved NDC such that it will be reimbursed by insurance. It is for the purpose of testing. By including the dilutant section **1802** in the testing kit, the NDC will still be associated with a “testing” kit and, therefore, it will be reimbursed as such. Even though it may be utilized to fill the prescription, it is still a “testing” kit. Additionally, each kit, regardless of what the purpose of it is, can have an NDC associated therewith and still be a testing kit. Thus, each regional kit will have NDC, each standardized kit will have NDC, and each combines standardized/regional kit will have NDC. In general, the term “regional kit” will refer to the kit that provides both standardized and regional antigens. The regional-only kit referred to as “regional-only.”

**[0093]** A group of regional kits can be provided, as illustrated in FIG. 22, one for each region that has defined issues. FIG. 23 illustrates a map of the United States and various regions that have differing requirements. They are the Pacific region, the Desert region, the Mountain region, Southwest region, the Plains region, the Great Lakes region, the South region and the Northeast region of the country. In the Pacific region, for example, trees usually pollinate from February to June. Trees that produce pollen that may trigger nasal allergy symptoms are cedar, walnut, and rye. Other trees that may cause nasal allergy symptoms in this region are Hazel, Juniper, Alder, Ash, Birch, Box Elder, and Oak. Grasses generally pollinate from March to November in this region. Some of the most common grasses that can trigger seasonal nasal allergy symptoms are sweet vernal, fescue, bluegrass, Bermuda grass, brome, orchard, quack grass, wheat grass, and redtop. The weed season for the region generally lasts from April to November, including ragweed, chenopod, and sage. Other weeds that can trigger nasal allergy symptoms include pigweed, iodine bush, saltbush, lamb’s quarter, Mexican fire bush, and Russian thistle.

**[0094]** In the Mountain region, trees usually pollinate from March to May. The trees that can trigger nasal allergy symptoms include the cedar, maple, and oak. Other trees that may trigger nasal allergy symptoms are the box elder, alder, birch, juniper, oak, and ash. Grasses pollinate from April to July in this region. Some of the grasses that can cause nasal allergy symptoms include quack grass/wheat grass, redtop, brome, Bermuda grass, and orchard grass. Weeds take over from early June through October, including ragweed, tumbleweed, and chenopod. Other weeds that are allergens or trigger nasal allergy symptoms include waterhemp, pigweed, iodine bush, saltbush/scale, sugarbeet, lamb’s quarter, Mexican fire bush, and Russian thistle. In the Desert region, cedar, ash, and oak pollinate from February to April. Other trees that can trigger nasal allergy symptoms include cypress, mesquite, ash, and olive. The grass season for this region can last from April until October. Some of the grasses that can cause nasal allergy symptoms include brome, Bermuda, salt grass, rye, Canary grass, and June grass. Weeds in the region pollinate from March to December, including ragweed, chenopod, and sage. Other weeds that are allergens or trigger nasal allergy symptoms include carelessweed, iodine bush, saltbush/scale, lamb’s quarter, and Russian thistle. In the Plains region, trees pollinate from March to May, including oak and cedar. Other trees that can trigger nasal allergy symptoms include alder, birch, maple/box elder, hazelnut, and oak. Grasses pollinate from late May to July in this region, including quack grass/wheat grass, redtop, brome, orchard grass, and rye grass. Weeds take over from July through September, including ragweed, plantain, and nettle. Other weeds that can trigger nasal allergy symptoms include water hemp, pigweed, lamb’s quarter, Mexican fire bush, and Russian thistle. In the Southwestern region, cedar, elm, and oak begin pollinating in early January and end in June. Other trees whose pollen is a known allergen include box elder and mesquite. Elm trees begin pollinating in late January and end in April, while olive pollinates from mid-March to mid-April.

**[0095]** The grass season for this region usually begins in April and lasts until September. The primary grasses for the region that can trigger nasal allergy symptoms are Bermuda, quack grass/wheat grass, and redtop. Weeds in the region include ragweed, chenopod, and dock, which pollinate

between July and October. Other weeds that can trigger nasal allergy symptoms include waterhemp, carelessweed/pigweed, saltbush/scale, and lamb's quarter. In the Great Lakes region, trees with pollen that can trigger nasal allergy symptoms are elm, maple/box elder, alder, birch, oak, and hickory. They can pollinate from March through June. Oak trees may continue pollinating until mid-June. Many grasses pollinate in this region during the summer, from May to July. Starting in May, redtop, brome, orchard grass, fescue, rye, Bermuda, and timothy pollinate. Weeds generally pollinate from July to September, including ragweed, plantain, and nettle. Other weeds that are allergens or trigger nasal allergy symptoms include waterhemp, lamb's quarter, pigweed, Mexican fire bush, and Russian thistle. Later in the season, hemp can trigger symptoms for people with seasonal nasal allergies—typically from mid-July to mid-August.

[0096] Then, from early August until late September, ragweed pollinates in this region. In this region, a variety of trees that can produce pollen that can trigger nasal allergy symptoms are pecan, oak, and cedar. Others include maple/box elder, birch, hickory, and oak. Grass pollens can pollinate nearly year-round in parts of this region. Some grasses that can trigger nasal allergy symptoms are Bermuda, saltgrass, bahia, redtop, vernal, orchard, rye, salt grass, fescue, and timothy. Among the weeds whose pollen can cause nasal allergy symptoms are ragweed, plantain, and nettle. Other weeds that are allergens or trigger nasal allergy symptoms are the lamb's quarter, sagebrush, English plantain, pigweed, waterhemp, and carelessweed. In the Northeast region, beginning in February and lasting until June, several types of trees in this region produce pollen that can trigger nasal allergy symptoms—particularly birch, maple/box elder, oak, juniper/cedar, and pine trees. From May to late August, grasses pollinate in the area including orchard, redtop, fescue, and timothy. Other grasses in the region that can trigger nasal allergy symptoms include vernal grass and Bermuda grass. From August to October, weeds such as ragweed, plantain, and nettle pollinate in this region. Other weeds that are allergens or trigger nasal allergy symptoms include lamb's quarter, cocklebur, pigweed, and Mexican fire bush.

[0097] It can be seen, therefore, that each region can have different requirements and each of these regional kits can be designed for a particular region and given as separate NDC code that can be reimbursable.

[0098] FIG. 24 is a diagrammatic view illustrating the fill operation of a sterile, single use vial of allergen with fill apparatus 2400. An allergen extract or concentrate is received from a vendor, usually in 5 ml bottles. The extract may, of course, be provided in larger or smaller volumes depending upon the particular allergen, the supplier and other factors. Typically, a vendor will provide an extract for a single allergen in a single container. The extracts are normally prepared as 1:10 to 1:20 w/v (weight to volume ratio) concentrates in glycerin. The extract 2402 may be placed in a container or the extract may be drawn directly from the supply bottle 2404. A tube 2406 is connected to bottle 2404 to draw the extract to a metering device 2408. In some variations, supply bottle 2404 may be provided with an injection stopper (not shown) and a needle is connected to tube 2406 and inserted into bottle 2404. A predetermined volume of extract is dispensed into sterile vials 2410 with a metering device 2408.

[0099] Metering device 2408 may be a small volume positive displacement pump, precision metering valve or similar device capable of repeatedly and accurately dispensing small volumes of extract 2402 into vials 2410. Vials 2410 may be glass or plastic and generally have a small volume, 1 ml or less. Prior to the filling operation, apparatus 2400 and vials 2410 are sterilized or provided in a sterile condition.

[0100] Referring to FIG. 25, after the extract has been dispensed into vials 2410, a sterile injection type stopper 2412 may be used to seal the vials. Stopper 2412 may fit into the top of the neck 2416 of the vial or simply extend over the mouth or opening 2418 of the vial. Since vials 2410 are designed for a single use, stoppers 2412 need not be as robust as the injection stoppers used for other applications where the stopper is penetrated multiple times with a syringe. Consequently stoppers 2412 may be a relatively thin layer of an appropriate material such as rubber, a plastic film or similar material that may tear or otherwise degrade when penetrated with a test prick or lancet. Stopper 2412 may fit into the top of the neck of the vial or simply extend over the mouth of the vial. In one embodiment, stoppers 2412 may be color-coded to indicate the specific allergen contained in vials 2410.

[0101] Turning to FIG. 26, after vials 2410 have been sealed with injection type stoppers 2412, a plastic or metal seal 2420 may be used to further secure the vial. For example, an aluminum peel or tear off seal 2420 may be crimped over injection type stopper 2412 and the mouth of vial 2410. Alternatively, a plastic snap off seal may be crimped over the injection type stopper 2412 to preserve the integrity of vial 2410. A combination seal, with an injection type stopper 2412 and an aluminum or plastic seal 2420 may be applied to vials 2410 in a single step. The filling and sealing operations are conducted under sterile conditions with sterile vials, stoppers and seals to minimize the possibility of contamination of the allergen extract. Seals 2420 may be color-coded, embossed or printed to indicate the specific allergen contained in vials 2410.

[0102] Turning to FIG. 27, after vials 2410 have been filled and sealed, vials 2410 may then be labeled. As illustrated, a label 2422 may be manually or mechanically applied to body of vial 2410. Label 2422 may include the specific allergen, the concentration of the allergen, the expiration date and other relevant information.

[0103] Turning to FIG. 28, after labeling, vials 2410 are packaged. In one embodiment, vials 2410 are mechanically or manually inserted into a plastic strip or tray 2424. Vials 2410 packaged in tray 2424 may then be boxed or placed in an appropriate container for transport. The vials 2410 may be used in a prick test by using a prick tester having a plurality of needles that are inserted through the stoppers 2412 of the vials 2410. The needles would then be applied to the skin of a patient to elicit a response. The prick tester may be placed on the vials 2410 as they reside in the tray 2424, or the vials 2410 may be placed in a different holding device having recesses for the vials, so that the body of the vials 2410 sit within the recesses and the stoppers 2412 protrude out from the recesses, similar to the prick tester shown in FIG. 10. Given that FIG. 10 illustrates a prick tester having wells instead of recesses for holding bottles, other prick tester configurations may be used such as those shown in U.S. patent application Ser. No. 15/222,709, published as U.S. Patent Application Publication No. 2017/

0027494, which U.S. patent application Ser. No. 15/222,709 and U.S. Patent Application Publication No. 2017/0027494 are herein incorporated by reference in their entirety. This prick testing kits hold vials in recesses arranged so that needles on a surface of the prick tester may be inserted into each vial and then applied to the skin of a patient. The prick testing kit may be assembled as one complete product, with the tray for the vials and the prick tester having the plurality of needles packaged together.

**[0104]** FIG. 29 is a flowchart illustrating one method of producing a sterile, single use vial of allergen for use in skin prick testing. The processes begins at **2900** when a source container of allergen extract, for example a 5 ml container of a 1:10 to 1:20 w/v (weight to volume ratio) allergen in glycerin is connected to a metering device such as a small volume positive displacement pump or similar metering device. At **2902** the metering device is used to dispense a small, for example 0.125 ml, volume of the allergen extract or concentrate into a sterile vial, for example, a 1 or 2 ml glass or plastic vial. An injection-type stopper is installed on the vial at step **2904** and a seal, for example, an aluminum or plastic peel off seal is crimped or otherwise applied over the injection-type stopper at **2906** to further secure the vial. In one embodiment, the filling, stoppering and sealing operations are performed with automated filling and sealing machines under sterile conditions to minimize the possibility of contamination. A label is applied to the vial at step **2908** and the filled, stoppered, sealed and labeled vial is packaged at step **2910**.

**[0105]** FIG. 30 is partial cross-sectional view of a filed, stoppered and labeled vial **2410**. A prick test normally requires only a small volume of extract, usually just enough to cover the tip of the lancet or needle used for the test. Consequently, the volume of extract required in the single use vials **2410** is small. Since the volume of extract required is small, vials **2410** are preferably filled with only enough extract to coat the lancet or needle for the test. For example, 0.1 or 0.2 ml of allergen extract may be dispensed into vial **2410**. The labeled, sterile single use vial **2410** of allergen extract may be packaged in a tray, strip, or boxed for distribution.

**[0106]** Although the preferred embodiment has been described in detail, it should be understood that various changes, substitutions and alterations can be made therein without departing from the spirit and scope of the invention as defined by the appended claims.

What is claimed is:

1. A method for producing prick test single use sterile vials for use in a prick test, comprising:

- providing at least one allergen source including there-within allergen extract;
- providing a plurality of single dose vials;
- connecting a tube to the at least one allergen source;
- connecting a metering device to the tube;
- drawing by the tube an amount of the allergen extract from one of the at least one allergen source;
- receiving the amount of allergen extract at the metering device; and
- dispensing by the metering device a volume of allergen extract into one of the plurality of single dose vials.

2. The method of claim 1, further comprising:  
sealing the one of the plurality of single dose vials with one of a plurality of stoppers; and

securing the one of the plurality of stoppers on the one of the plurality of single dose vials with one of a plurality of seals.

3. The method of claim 2, further comprising affixing one of a plurality of labels to the one of the plurality of single dose vials.

4. The method of claim 3, further comprising:  
providing a tray having disposed therein a plurality of recesses; and

placing within one of the plurality of recesses of the tray the one of the plurality of single dose vials.

5. The method of claim 4, further comprising repeating the drawing, receiving, dispensing, sealing, securing, affixing, and placing steps until a desired number of the plurality of single dose vials are produced.

6. The method of claim 5, further comprising:  
providing a prick tester having a plurality of needles extending thereon;

aligning the plurality of needles of the prick tester with the plurality of single dose vials;

inserting each of the plurality of needles of the prick tester into the plurality of single dose vials; and

applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

7. The method of claim 5, further comprising:  
extracting a predetermined amount of allergen extract from one of the plurality of single dose vials; and  
disposing the predetermined amount of allergen extract into a corresponding one of a plurality of wells.

8. The method of claim 7, further comprising repeating the extracting and disposing steps until a desired number of the plurality of wells contain allergen extract.

9. The method of claim 8, further comprising:  
providing a prick tester having a plurality of needles extending thereon;

aligning the plurality of needles of the prick tester with the plurality of wells;

inserting each of the plurality of needles of the prick tester into one of the plurality of wells; and

applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

10. The method of claim 1, further comprising:  
providing a prick tester having a plurality of needles extending thereon;

aligning the plurality of needles of the prick tester with the plurality of single dose vials;

inserting each of the plurality of needles of the prick tester into the plurality of single dose vials; and

applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

11. The method of claim 1, further comprising:  
extracting a predetermined amount of allergen extract from one of the plurality of single dose vials; and  
disposing the predetermined amount of allergen extract into a corresponding one of a plurality of wells.

12. The method of claim 11, further comprising repeating the extracting and disposing steps until a desired number of the plurality of wells contain allergen extract.

13. The method of claim 12, further comprising:  
providing a prick tester having a plurality of needles extending thereon;

aligning the plurality of needles of the prick tester with the plurality of wells;

inserting each of the plurality of needles of the prick tester into one of the plurality of wells; and  
applying the plurality of needles of the prick tester to the skin of a patient to elicit a potential response.

**14.** A system for producing prick test single use sterile vials for use in a prick test, comprising:

- an allergen source including therewithin allergen extract;
- a tube connected to the allergen source to draw the allergen extract from the allergen source;
- a metering device connected to the tube, wherein the metering device receives by way of the tube an amount of allergen extract from the allergen extract in allergen source, and wherein the metering device serves to dispense a volume of extract;
- a plurality of single dose vials, wherein each one of the plurality of single dose vials is configured to receive the volume of extract dispensed by the metering device;
- a plurality of sterile injection stoppers for use in sealing the plurality of single dose vials;
- a plurality of seals to secure the plurality of sterile injection stoppers on the plurality of single dose vials;
- a plurality of labels, each one of the plurality of labels configured to be affixed to one of the plurality of single dose vials;

a tray having disposed therein a plurality of recesses, each one of the plurality of recesses able to receive one of the plurality of single dose vials; and

a prick tester having a plurality of needles extending thereon, wherein the plurality of needles are arranged to be inserted into the plurality of single dose vials disposed in the plurality of recesses of the tray.

**15.** The system of claim **14**, further comprising a dilutant for dispensing into a number of the plurality of single dose vials.

**16.** The system of claim **15**, wherein the dilutant comprises a saline solution.

**17.** The system of claim **15**, wherein the dilutant comprises a glycerol solution.

**18.** The system of claim **14**, further comprising a first control material and a second control material, wherein the first control material is dispensed into a number of the plurality of single dose vials, and wherein the second control material is dispensed into a number of the plurality of single dose vials.

**19.** The system of claim **18**, wherein the first control material is water.

**20.** The system of claim **19**, wherein the second control material is a histamine.

\* \* \* \* \*